


Welcome Guest from United States

[Sign In](#) | [Change Country](#)
 0 Items
Search: 
[Product Browse](#) [Products >>... >> Meters >> Benchtop pH meters](#)

## Orion 920Aplus Multi-Channel Benchtop Meter

- Advanced Features for ISE Measurements
- Automatic Blank Correction Calculations
- Electrode Switchbox compatible for Up to 7 channels
- Three (3) year meter warranty

### Product Detail

The Orion 920Aplus provides more features than ever before, making this meter easy to use and smart. Meter accepts up to five standards for concentration calibration to ensure low-level measurement accuracy. Autoblank correction feature automatically corrects concentration results for background ion levels. Attach the Orion 607 Electrode Switchbox to work with up to seven different electrodes.

Other features adding to the ease of use include high and low alarm limits and display hold. Orion 920Aplus advanced ISE features include pre-programmed incremental techniques, calibration edit, and temperature-compensated ISE measurement. Use calibration edit feature for reviewing calibration standards to determine working range. With GLPdoc™, all measurement parameters, including slope, temperature, mV of standards, and time and date of calibration, automatically print on the optional Orion 900A printer or to a computer via RS232.

### Orion 920Aplus Multi-Channel Benchtop Meter


Product#	Price	Package	
<b>0920A0</b>	-	110 V with pH Electrode	<a href="#">SELECT</a>
<b>0920A3</b>	-	110 V with Electrode Holder/Stand	<a href="#">SELECT</a>
<b>0920A9</b>	-	110 V Meter and Adapter Only	<a href="#">SELECT</a>
<b>0920A1</b>	-	220 V with pH Electrode	<a href="#">SELECT</a>
<b>0920A8</b>	-	220 V Meter and Adapter only	<a href="#">SELECT</a>
<b>0920A7</b>	-	240 V Meter and Adapter only	<a href="#">SELECT</a>
<b>0920A2</b>	-	Meter and Manual Only	<a href="#">SELECT</a>

[Details All Product Numbers](#)
[Details All Product Numbers](#)

### Purchase Details

 [Request Quote](#)

### Product Contact

 [Sales Contact](#) |  [Service Contact](#)
 [North America +1 800 225 1480](#)

### Literature Request

### Postal Mail:

- Orion Laboratory Products Catalog

### Specifications:

<b>pH Range :</b>	-2.000 to 19.999
<b>pH Resolution :</b>	0.001/0.01/0.1
<b>pH Relative Accuracy :</b>	± 0.002
<b>pH Slope :</b>	80 to 120%
<b>Auto-buffer Recognition :</b>	1.68, 4.01, 7.00, 10.01, 12.46
<b>Calibration Points :</b>	5
<b>Concentration Range :</b>	0.000 to 19900
<b>Concentration Resolution :</b>	± 1 least significant digit
<b>Concentration Relative Accuracy :</b>	± 0.5% of reading
<b>Temperature Range :</b>	-5.0 to 105.0°C
<b>Temperature Resolution :</b>	0.1°C
<b>Temperature Relative Accuracy :</b>	± 0.1°C
<b>Millivolt Range :</b>	± 1600.0 mV
<b>Millivolt Resolution :</b>	0.1 mV
<b>Millivolt Relative Accuracy :</b>	± 0.2 mV or ± 0.05% of reading, whichever is greater
<b>Relative Millivolt Range :</b>	± 1999.9 mV
<b>Relative Millivolt Resolution :</b>	0.1 mV
<b>Display :</b>	Custom color vacuum fluorescent
<b>Inputs :</b>	2 BNC, 2 pin tip, ATC, power, RS232
<b>Outputs :</b>	RS232, 2.5 mm recorder, Karl Fischer
<b>Power Requirements :</b>	AC Line 110V, 115V, 220V or 240 V
<b>Agency Approvals :</b>	CE, CSA, UL
<b>Environmental Requirements :</b>	5°C to 45°C and 5 to 85% relative humidity, non-condensing