

A range of high quality and precision burettes at a reasonable price

They cover the applications between the manual digital dispensers and the expensive motorized burettes.

Possibility to be supplied with 1 or 2 additional peristaltic pumps.

The piston's stroke is divided in 40000 steps.

The resolution on display is 1 µl.

The main applications are:

- Solution preparation.
- Repetitive reagent dispensing.
- Standard dilution.
- Titrations with optical indicator.
- etc.



Model	Burette 15 DPD	Burette 1S	Multi-Burette 2S	Multi-Burette 4S
Number of syringes	1	1	2	4
Magnetic stirrer (optional)	1	1	2	4
Dilution programs	Yes	—	—	—
Dispensing programs	—	—	Yes	—
Preparation programs	—	—	—	Yes
Additional peristaltic pumps	—	—	optional, 1 or 2	—